

# Evaluating a Fully Automated 16S/tNGS Solution for the Detection of Bacterial and Fungal Pathogens in Formalin-fixed, Paraffin-embedded (FFPE) tissues


R. Yee<sup>1</sup>, J. H. J. Ng<sup>2</sup>, L. Codoy<sup>3</sup>, H. Simonyan<sup>1</sup>, O. Kochar<sup>3</sup>, N. Ordaz-Yuan<sup>2</sup>, A. Paropkari<sup>2</sup>, D. Tran<sup>2</sup>, A. Lin<sup>2</sup>, R. Khaksar<sup>2</sup>

<sup>1</sup> Department of Pathology, George Washington University School of Medicine and Health Sciences, Washington, DC, United States, <sup>2</sup> Clear Labs, Inc., San Carlos, CA, United States, <sup>3</sup>Department of Laboratory Medicine, George Washington University Hospital, Washington, DC, United States

## Introduction

Formalin-fixed, paraffin-embedded (FFPE) tissue specimens in pathology can preserve the integrity of the tissue and allow for histopathological examination to aid in the diagnosis of infectious diseases. However, exact determination of the infectious organism can be difficult. Surgical pathology may reveal morphological features suggestive of infection, but these findings alone may be inconclusive due to nonspecific staining or morphological similarities among organisms. Next generation sequencing (NGS), specifically targeted 16S rRNA and ITS gene sequencing (tNGS) has become a widely used approach for identification of various bacteria and fungi, respectively. A streamlined, automated tNGS workflow has the potential to overcome these challenges, providing higher sensitivity and broader pathogen coverage while minimizing manual steps and turnaround time. Here, we evaluated the performance of a commercially-available, fully automated 16S/ITS tNGS test on FFPE tissues of various anatomic sites.

## Methods




### Assessment:

- Determination of concordance, sensitivity, specificity, positive and negative predictive values using SoC results as gold standard

## Results

### A: Anatomic Sources of FFPE Tested

A total of 30 specimens were tested from various anatomic sites.



### B: Performance of Clear Dx tNGS

	Overall	Bacterial	Fungal
Sensitivity	10/12 (83.3%)	6/8 (75%)	4/4 (100%)
Specificity	18/18 (100%)	22/22 (100%)	26/26 (100%)
Positive Predictive Value	10/10 (100%)	6/6 (100%)	4/4 (100%)
Negative Predictive Value	18/20 (90%)	22/24 (92%)	26/26 (100%)


### C: Comparison of Results Between tNGS and Pathology Review

Histopathology Observation	tNGS Concordance
<b>Positive (total n= 12)</b>	10/12 (83.3%)
<b>Bacterial (n=8)</b>	6/8 (75.0%)
• Positive for bacterial organisms	2/2 (100%)
• Positive for intestinal spirochetosis	4/4 (100%)
• Positive for mycobacterial organisms	0/1 (0%)
• Positive for <i>H. pylori</i>	0/1 (0%)
<b>Fungal (n=4)</b>	4/4 (100%)
• Positive for fungal organism (e.g., yeast, spores)	4/4 (100%)
<b>Negative (n=18)</b>	18/18 (100%)
<b>Grand total (N=30)</b>	28/30 (93.3%)

- Clear Dx™ Microbial ID tNGS achieved an overall concordance of 93.3%.
- 100% concordance for detection of yeast fungus.
- 75% concordance for bacterial organisms.
- tNGS missed two cases where SoC reported positive for acid-fast bacilli and *H. pylori*.

NGS Target	Read Abundance (raw) [mean, (range)]	Read Abundance (%) [mean, (range)]
<b>Overall</b>	70,205 (2,141 – 206,460)	32.38 (4.4 – 88.6)
<b>Bacteria</b>		
• Vaginal Flora	29,769 (13,757-73,487)	9.57 (4.4-23.6)
• <i>Pseudomonas</i> species	17,736	22.96
• <i>Brachyspira</i> species	81,193 (2,141-144,542)	71.08 (42.5 – 88.6)
<b>Fungal</b>		
• <i>Candida albicans</i>	113,133 (46,218 – 206,460)	24.12 (6.0-48.9)
• <i>Nakaseomyces</i> species	10,367	14.34


### D: Examples of Histopathological Findings and Respective tNGS Results



GMS stain

#### Stomach Biopsy


- Path report: Positive for fungal organism consistent with *Candida*
- tNGS: *Candida albicans* (reads: 46218; 6%)



PAS stain

#### Stomach Biopsy


- Path report: Positive for fungal spores
- tNGS: *Nakaseomyces* spp. (reads: 10,367; 4%)



Gram stain

#### Cervix Biopsy

- Path report: Degenerative squamous epithelium with bacteria
- tNGS: Vaginal flora [*Lactobacillus* spp., (reads: 73,487; 24%); *Achromobacter* spp. (reads: 15,573; 5%); *Cutibacterium* spp. (reads: 13,757; 4%)]



Warthin-Starry

#### Stomach Biopsy

- Path report: Positive for intestinal spirochetosis, negative for *H. pylori*
- tNGS: *Brachyspira aalborgi* (reads: 96,879; 82%)

## Conclusions

- The fully automated Clear Dx™ Microbial ID tNGS assay is a promising approach for accurate and unbiased fungal and bacterial identification from FFPE
- Performance of Clear Dx™ Microbial ID tNGS assay is comparable from various anatomical sites.
- Ongoing runs with additional FFPE specimens will increase the power of this study.
- Further evaluation of extraction protocols and inclusion of specimens with more diverse pathological findings (e.g., fungal hyphae) will provide insights into the accuracy.
- Sequencing approaches with appropriate diagnostic stewardship can complement SOC practices and will have a significant impact in infectious disease management decisions.

### References:

- Larkin *et al.*, Amplicon-Based Next-Generation Sequencing for Detection of Fungi in Formalin-Fixed, Paraffin-Embedded Tissues: Correlation with Histopathology and Clinical Applications. *J Mol Diagn.* 2020 Oct;22(10):1287-1293. doi: 10.1016/j.jmoldx.2020.06.017. Epub 2020 Jul 30. PMID: 32738297.
- Lao H-Y, et al. The clinical utility of Nanopore 16S rRNA gene sequencing for direct bacterial identification in normally sterile body fluids. *Frontiers in Microbiology.* 2024;14:1324494. DOI: 10.3389/fmicb.2023.1324494.
- Flurin L, et al. Targeted Metagenomic Sequencing-based Approach Applied to 2146 Tissue and Body Fluid Samples in Routine Clinical Practice. *Clin Infect Dis.* 2022 Nov 14;75(10):1800-1808. doi: 10.1093/cid/ciac247. PMID: 35362534; PMCID: PMC9662179.