# Clear Labs Clear Safety<sup>TM</sup> Salmonella: Automated Targeted NGS Detection and Serotyping from Sample Enrichment

Andrew Lin, Atul K. Singh, Ramin Khaksar

# BACKGROUND

Salmonella serotyping is an important tool to help identify strains of concern for implementation of risk-based mitigation strategies. Current methods (traditional anti-sera, DNA arrays, bead-based arrays, and whole genome sequencing) require that analysts culture and isolate Salmonella colonies prior to performing the serotyping methods.

Clear Safety<sup>™</sup> Salmonella uses an automated NGS-based platform to simultaneously detect Salmonella and identify the most common serotypes from sample enrichments avoiding the need for isolating Salmonella colonies and reducing the time to results from >72h to 40h.

Day 1	Day 2	Day 3	Day 4
6h12h18h24hSample Enrichment	30h 36h 42h 48h Secondary Enrichment	54h 60h 66h 72h Selective Agars	78h 84h 90h 96h Serotyping/ Confirmation
Sample Enrichment	Clear Safety		

# Highlighted Features of Clear Safety Salmonella

- Targeted NGS generates millions of sequencing reads to detect multiple genes for each target pathogen that allows for built-in redundancy resulting in less false negatives and false positives.
- Automated platform simplifies the workflow minimizing user intervention in performing the assay and result analysis.
- NGS technology allows for high throughput (192) samples) detection of Salmonella and serotyping of the 63 most common Salmonella serotypes
- Clear Safety Salmonella targets DNA from live cells, reducing false positives from dead cell DNA
- Clear Safety Salmonella is AOAC and NPIP-approved for detection in food and environmental samples


## Figure 1: Clear Safety Salmonella: An Automated, **Targeted NGS Platform**



Integrated, end-to-end automation platform for direct detection and serotyping of Salmonella from sample enrichment. Steps include sample lysis, live and dead sample treatment, PCR, library preparation, sequencing and analytics



For more information regarding Clear Safety Salmonella, please contact us at inquiries@clearlabs.com



# Sampling

A total of 363 poultry environmental samples (boot swabs, dust swabs and hatchery waste) sampled from primary production houses were collected following the National Poultry Improvement Plan (NPIP) program standards. Swabs were either pre-moistened or had 30mL of double strength skim milk added to them. Samples were analyzed by three independent laboratories following NPIP recommendations. Samples were stored at refrigerator temperatures (2-4°C) for no more than 5 days.

# Enrichment

Poultry environmental samples were enriched in 150mL (1:10 sample to media ratio) of pre-warmed tetrathionate broth (TT) and incubated at  $37 \pm 2^{\circ}$ C for 28-48h.

### Clear Safety Salmonella

- 1. 50µL of TT enrichment was added to 450µL CL Prep Solution in a 1.2 mL sample tube in a sample rack.
- 2. Sample rack and reagents were loaded onto the Clear Safety Platform. Bacteria lysis, PCR, library preparation and sequencing were performed by the automated workflow.

### **NPIP** Reference Method

- 1. 100µL of TT enrichment were spotted on MSRV incubated at  $42 \pm 2^{\circ}$ C for 24h
- 2. Outer edge of MSRV growth were streaked on BGN and XLT4 agar, incubated at 35 ± 2°C for 24h
- 3. Suspect colonies were streaked on TSI and LIA slants incubated at  $35 \pm 2^{\circ}$ C for 24h.
- 4. Salmonella isolates were serotyped with anti-sera.

In this paired study, 363 poultry environmental samples were analyzed by Clear Safety Salmonella and the NPIP program standard. With a Cohen's Kappa measure of inter-rater reliability of 0.95, there is **near perfect agreement** of the two methods for detection of Salmonella. Clear Safety Salmonella proves to be accurate with diagnostic sensitivity, specificity, positive predictive value and negative predictive values all surpassing the USDA MLG acceptance criteria for alternative methods of > than 90%.

Di Di

Clear Safety Salmonella serotyping results from sample enrichments agreed with traditional anti-sera serotyping from isolates in 86.1% of the samples.

Many factors may have contributed to the lack of agreement in serotyping results for some samples. While classical anti-sera serotyping is included in the NPIP reference method, the accuracy of this method can vary depending on the quality of anti-sera and experience of analysts. In addition, Clear Safety Salmonella results from enrichment could differ from classical serotyping of isolates if the sample contained multiple serotypes and the analyst did not pick enough isolates to get a full representation. Bias in secondary selective enrichments favoring growth of certain serotypes over others could also contribute to these discrepancies.

# RESULTS

Table 1: Consolidated	Data from	<b>Coordinating Labs</b>
-----------------------	-----------	--------------------------

<b>NPIP Reference Method</b>		Totolo
Positive	Negative	Totals
113	3	116
5	242	247
118	245	363
	Positive 113 5	PositiveNegative11335242

Statistical parameter	Consolidated data	Acceptance criteria
Cohen Kappa	0.95	0.81 - 1.00
iagnostic sensitivity (dSN)	95.8	> 90%
iagnostic specificity (dSP)	98.8	> 90%
Positive Predictive Value (PPV)	95.8	> 90%
Vegative Predictive Value (NPV)	98.0	> 90%

# Table 2: Statistical Analysis of the Data

### Table 3: Clear Safety Serotyping vs Classical Anti sera

Laborat Laborat Laborat Tota

# CONCLUSION

Clear Safety Salmonella is an automated sequencing system that accurately detects Salmonella and confirms serotypes 12h after enrichment, compared to the culture confirmation methods that takes around 4-7 days. In this independent study, the Clear Safety Salmonella assay has shown high agreement with the reference method. Thus, the findings of this study support the claim that the Clear Safety Salmonella assay presents a valuable, affordable, and high throughput automated solution that can improve identification of Salmonella in poultry production houses. Clear Safety Salmonella offers rapid detection and serotyping that helps food producers make informed decisions in implementing risk-based mitigation strategies to minimize Salmonella contamination and outbreaks.

# Table 4: List of Identifiable serotypes from Clear Safety Salmonella

Enteritidis Typhimur I 4,[5],12: Newport Kentucky Abaetetu Agona Alachua Albany Anatum Bareilly Berta Blockley

The Clear Safety Salmonella platform offers molecular sequence-based serotyping of 63 serotypes that is **98.4%** of most commonly identified in the poultry industry.

 $\Lambda$ 

© 2022 Clear Labs Inc.

1559 Industrial Road San Carlos, CA 94070

	No. of Positive Samples	Serotype Matched with Ref Method	Serotyping Accuracy
tory-1	48	40	83.3%
tory-2	35	30	85.7%
tory-3	25	23	92.0%
al	108	93	86.1%

is	Braenderup	Muenchen	Virchow	Gaminara
irium	Cerro	Oranienburg	Tennessee	Havana
2:i:-	Derby	Panama	Liverpool	Idikan
	Dublin	Paratyphi B	Uganda	Lille
У	Give	Poona	Gallinarum/Pullorum	Pomona
uba	Hadar	Reading	Minnesota	Putten
	Heidelberg	Rissen	Ohio	Roodepoort
	Infantis	Saintpaul	Molade	Muenster
	Javiana	Schwarzengrund	Litchfield	Norwich
	Johannesburg	Senftenberg	Sandiego	Worthington
	Mbandaka	Stanley	Meleagridis	Ouakam
	Mississippi	Thompson	Barranquilla	
,	Montevideo	Typhi	Cubana	